

MULTI CHANNEL STIMULATOR

MICRO-CONTROLLER BASED
ABS SHOCKPROOF BODY
AESTHETIC LOOK
COMPACT AND EASY TO USE
LIGHT WEIGHT
MEMBRANE KEY BOARD
4 LINE LCD DISPLAY
SIX TYPES OF CURRENTS
DIAGNOSTIC/THERAPETIC

TECHNICAL SPECIFICATIONS :

- 1) INTENSITY : 0 – 100, ADJUSTABLE
- 2) DISPLAY : 4 LINE LCD DISPLAY
- 3) MODES : 2 POLE AM, INTERRUPTED GALVANIC, FARADIC, SURGE FARADIC, TNS, ACUPUNCTURE
- 4) INT. GAL : FREQ : 1, 5, 10, 50, 100 Hz, DURATION : 0.01, 0.03, 0.1, 0.3, 1, 3, 10, 30, 100
- 5) SUR. FAR : ON TIME : 1 – 10 SEC, OFF TIME : 1 – 10 SEC
- 6) TNS FREQ : 2, 50, 80, 100, 150 Hz
- 7) 2 POLE AM : 5 KHz
- 8) ACUPUNCTURE : 2, 50, 80, 100, 150 Hz
- 9) FARADIC : 40 Hz
- 10) RESOLUTION : 1 VOLT
- 11) INPUT : 230VAC (MAINS OPERATED)
- 12) One Pair of Cable with Electrodes, Straps – 1 Pair
Handle with cable with electrode – 1 no. Carrying case – 1no
- 13) WEIGHT : < 2 Kg
- 14) DIMENSION : 275 x 200 x 85 mm



FOUR CHANNEL

MICRO-CONTROLLER BASED
AUTO CUT-OFF AFTER PRESET TIME

COMPACT EASY TO USE
20 PROGRAMS

AESTHETIC LOOK
SMOOTH OUTPUT CONTROL

TECHNICAL SPECIFICATIONS :

- 1) INTENSITY : 0 – 120 VOLTS, ADJUSTABLE
- 2) DISPLAY : 2 SEVEN SEGMENT DISPLAY
- 3) MODES : INT. GALVANIC, FARADIC, SUR. FARADIC
- 4) PROGRAMS : TEN FOR INT. GALVANIC
TEN FOR FARADIC
- 5) TIMER : 0 – 99MIN
- 6) INPUT : 230VAC (MAINS OPERATED - ADAPTOR)
- 7) One Pair of Cable with Electrodes, Straps – 1 Pair
- 8) CARRYING CASE
- 9) DIMENSION : 170 x 130 x 50 mm
- 10) WEIGHT : < 1 Kg



MINI MS - DOUBLE CHANNEL

GALTRON ELECTROMEDICAL PVT. LTD.

AN ISO 9001 : 2008 CERTIFIED ORGANISATION

G/I , Nahar & Seth Industrial Estate, Chakala Road, Andheri – E, Mumbai – 400099.

PH : 022 – 28345018, 28323118, 66910067. Email ID : galtron@rediffmail.com

Website : www.galtronindia.com

