

# MICRO-CONTROLLER BASED DIAGNOSTIC STIMULATOR WITH TNS

MICRO-CONTROLLER BASED  
COMPACT AND EASY TO USE  
2 LINE LCD DISPLAY

ABS SHOCKPROOF BODY  
LIGHT WEIGHT  
FIVE TYPES OF CURRENTS

AESTHETIC LOOK  
MEMBRANE KEY BOARD  
DIAGNOSTIC/THERAPETIC



MODEL : DSL - 007

## TECHNICAL SPECIFICATIONS :

- 1] INTENSITY : 0 – 120 VOLTS, ADJUSTABLE
- 2] DISPLAY : 2 LINE LCD DISPLAY
- 3] MODES : GALVANIC, INTERRUPTED GALVANIC, FARADIC, SURGE FARADIC, TNS
- 4] INT. GAL : FREQ : 1, 5, 10, 50, 100 Hz, DURATION : 0.01, 0.03, 0.1, 0.3, 1, 3, 10, 30, 100
- 5] SUR. FAR : ON TIME : 1 – 10 SEC, OFF TIME : 1 – 10 SEC
- 6] TNS FREQ : 2, 50, 80, 100, 150 Hz
- 7] RESOLUTION: 1 VOLT
- 8] INPUT : 230VAC ( MAINS OPERATED )
- 9] One Pair of Cable with Electrodes, Straps – 1 Pair  
Handle with cable with electrode – 1 no. Carrying case – 1no
- 10] WEIGHT : < 2 Kg
- 11] DIMENSION : 275 x 200 x 85 mm

# MINI MUSCLE STIMULATOR

MICRO-CONTROLLER BASED  
AUTO CUT-OFF AFTER PRESET TIME

COMPACT EASY TO USE  
20 PROGRAMS  
SMOOTH OUTPUT CONTROL



## TECHNICAL SPECIFICATIONS :

- 1] INTENSITY : 0 – 120 VOLTS, ADJUSTABLE
- 2] DISPLAY : 2 SEVEN SEGMENT DISPLAY
- 3] MODES : INT. GALVANIC, FARADIC, SUR. FARADIC
- 4] PROGRAMS : TEN FOR INT. GALVANIC  
TEN FOR FARADIC
- 5] TIMER : 0 – 99MIN
- 6] INPUT : 230VAC ( MAINS OPERATED )
- 7] One Pair of Cable with Electrodes, Straps – 1 Pair
- 8] CARRYING CASE
- 9] DIMENSION : 90 x 155 x 50 mm
- 10] WEIGHT : < 1 Kg

**Our Range of Products :** Electrosurgical cautery, (CTG) Fetal Monitors, Fetal Doppler's, Colposcope, Pulse Oximeters, Autoclaves, Suction, Radiant Heat Warmer, Phototherapy, Shortwave Diathermy, Traction, C.P.M, Muscle Stimulators, Interferential therapy unit, Ultrasonic Therapy unit, Skin Cautery, Cryo unit, Spotlight, etc.

**GALTRON ELECTROMEDICAL PVT. LTD.**

AN ISO 9001 : 2008 CERTIFIED ORGANISATION

G/I , Nahar & Seth Industrial Estate, Chakala Road, Andheri – E, Mumbai – 400099.

PH : 022 – 28345018, 28323118, 66910067. Email ID : [galtron@rediffmail.com](mailto:galtron@rediffmail.com)

Website : [www.galtronindia.com](http://www.galtronindia.com)

