



MICRO-CONTROLLER BASED NERVE STIMULATOR

MICRO-CONTROLLER BASED
COMPACT
MEMBRANE KEY BOARD
AESTHETIC LOOK

ABS SHOCKPROOF BODY
LIGHT WEIGHT
2 LINE LCD DISPLAY
EASY TO USE



TECHNICAL SPECIFICATIONS :

1) INTENSITY	:	0 – 120 , ADJUSTABLE
2) DISPLAY	:	2 LINE LCD DISPLAY
3) CURRENT	:	PULSE
4) PULSE FREQUENCY	:	1, 2, 5 Hz
		21 Programs
5) ELECTRODES	:	1 PAIR, Handle with cable with electrode – 1 no.
8) RESOLUTION	:	1 VOLT
9) INPUT	:	230VAC (MAINS OPERATED)
10) CARRYING CASE	:	1 NO.

Our Range of Products : Electrosurgical cautery, (CTG) Fetal Monitors, Fetal Doppler's, Colposcope, Pulse Oximeters, Autoclaves, Suction, Radiant Heat Warmer, Phototherapy, Shortwave Diathermy, Traction, C.P.M, Muscle Stimulators, Interferential therapy unit, Ultrasonic Therapy unit, Skin Cautery, Cryo unit, Spotlight, etc.

GALTRON ELECTROMEDICAL PVT. LTD.

AN ISO 9001 : 2008 CERTIFIED ORGANISATION

G/I , Nahar & Seth Industrial Estate, Chakala Road, Andheri – E, Mumbai – 400099.

PH : 022 – 28345018, 28323118, 66910067. Email ID : galtron@rediffmail.com

Website : www.galtronindia.com